

Product Datasheet

WDR5B Antibody (orb1270415)

Description

WDR5B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This WDR5B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 58-85 amino acids from the N-terminal region of human WDR5B.

Target

WDR5B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


36 kDa

Uniprot ID
[Q86VZ2](#)
NCBI
[Q86VZ2](#)
Dilution Range

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

WiDr
72


Western blot analysis in WiDr cell line

...



Flow cytometric analysis of WiDr cells (...)