

Product Datasheet

ICAM2 Antibody (orb1270420)

Description
ICAM2 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This ICAM2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 191-218 amino acids from the Central region of human ICAM2.

Target

ICAM2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


31 kDa

Uniprot ID
P13598
NCBI
P13598
Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis in HepG2 cell line...



Flow cytometric analysis of HepG2 cells

...