

## Product Datasheet

### RAB3IL1 Antibody (orb1270429)

**Description**

RAB3IL1 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

FC, WB

**Applications**
**Immunogen**

This RAB3IL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 297-325 amino acids from the C-terminal region of human RAB3IL1.

**Target**

RAB3IL1

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**

43 kDa


**Uniprot ID**
[Q8TBN0](#)
**NCBI**
[Q8TBN0](#)
**Dilution Range**

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50


**Expiration Date**

12 months from date of receipt.

NCI-H460



Western blot analysis in NCI-H460 cell l...



Flow cytometric analysis of NCI-H460 cel...