

Product Datasheet

TTC39B Antibody (orb1270441)

Description

TTC39B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This TTC39B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 625-654 amino acids from the C-terminal region of human TTC39B.

Target

TTC39B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


77 kDa

Uniprot ID
[Q5VTQ0](#)
NCBI
[Q5VTQ0](#)
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Western Blot at
1:1000 dilution + THP-
1 ...



Western blot analysis
in mouse testis ti...



TTC39B antibody
immunohistochemistry
ana...