



Product Datasheet

MCM2 antibody (orb127046)


Description	Rabbit Polyclonal to MCM2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P49736
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 2...



Western blot analysis of extracts from 3...



Immunohistochemical analysis of formalin...