

Product Datasheet

TNFAIP8 Antibody (orb1270473)

Description

TNFAIP8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications
Immunogen

This TNFAIP8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 50-78 amino acids from the N-terminal region of human TNFAIP8.

Target

TNFAIP8

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


23 kDa

Uniprot ID
[O95379](#)
NCBI
[O95379](#)
Dilution Range

For WB starting dilution is: 1:1000 For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.



Western blot analysis in K562 cell line

...



Flow cytometric analysis of K562 cells (...)