



Product Datasheet

TET1 antibody (orb12705)


Description	Rabbit polyclonal antibody to TET1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human TET1(NP_085128.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	250KDa
Uniprot ID	Q8NFU7
Dilution Range	WB: 1:500-1:2000, IP: 1:20-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of mouse te...



Immunohistochemical staining of mouse sp...