

## Product Datasheet

### CCDC126 Antibody (orb1270542)

## Description

CCDC126 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

This CCDC126 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 77-106 amino acids from the C-terminal region of human CCDC126.

### Target

CCDC126

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

16 kDa

### Uniprot ID

[Q96EE4](#)

### NCBI


[Q96EE4](#)

### Dilution Range

For WB starting dilution is: 1:1000

### Expiration Date

12 months from date of receipt.



Western blot analysis in mouse liver tis...