

Product Datasheet

ACAD10 Antibody (orb1270571)

Description

ACAD10 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This ACD10 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 393-421 amino acids from the Central region of human ACD10.

Target

ACAD10

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


119 kDa

Uniprot ID
Q6JQN1
NCBI
Q6JQN1
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

NCI-H460
250


Western blot analysis
in NCI-H460 cell l...



ACD10 Antibody
immunohistochemistry
anal...



Flow cytometric
analysis of NCI-H460
cel...