

Product Datasheet

RPH3AL Antibody (orb1270623)

Description

RPH3AL Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This RPH3AL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 208-237 amino acids from the C-terminal region of human RPH3AL.

Target

RPH3AL

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

[Q9UNE2](#)

NCBI


[Q9UNE2](#)

Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.



Western Blot at
1:1000 dilution +
A549 w...



RPH3AL Antibody
immunohistochemistry
ana...