

Product Datasheet

ITGAM Antibody (orb1270680)

Description

ITGAM Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This CD11b antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 253-282 amino acids from the N-terminal region of human CD11b.

Target

ITGAM

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:500For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


127 kDa

Uniprot ID
P11215
NCBI
P11215
Dilution Range


For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:500For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:2000
dilution +
human ...



Antibody
staining
D11b in
human tonsil
t...



Antibody
staining
D11b in
human
spleen t...