

Product Datasheet

ARV1 Antibody (orb1270688)

Description

ARV1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

FC, IHC-P, WB

Immunogen

This ARV1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 23-51 amino acids from the N-terminal region of human ARV1.

Target

ARV1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

31 kDa

Uniprot ID

[Q9H2C2](#)

NCBI


[Q9H2C2](#)

Dilution Range


For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Overlay histogram showing A2058 cells stained with ARV1 antibody.



Overlay histogram showing A2058 cells stained with ARV1 antibody.



Western Blot analysis of ARV1 expression in various cell lines. Lane 1: ...