

## Product Datasheet

### C1QTNF9B-AS1 Antibody (orb1270698)

## Description

C1QTNF9B-AS1 Antibody

## Species/Host

Rabbit

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

FC, IHC-P, WB

## Applications

## Immunogen

This PCOTH antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 14-42 amino acids from the Central region of human PCOTH.

## Target

C1QTNF9B-AS1

## Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

## Form/Appearance

Liquid

## Concentration

batch dependent

## Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

## Note

For research use only

## Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

## Isotype

Rabbit Ig

## Clonality

Polyclonal

## MW

11 kDa

## Uniprot ID

[Q58A44](#)

## NCBI


[Q58A44](#)

## Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

## Expiration Date


12 months from date of receipt.



Flow cytometric analysis of 293 cells (r...



Western blot analysis in 293 cell line l...



PCOTH antibody immunohistochemistry anal...