


Product Datasheet

TBCCD1 Antibody (orb1270736)

Description	TBCCD1 Antibody		
Species/Host	Rabbit		
Reactivity	Human, Mouse		
Conjugation	Unconjugated		
Tested Applications	WB	Western blot analysis in mouse liver tis...	
Immunogen	This TBCCD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 491-520 amino acids from the C-terminal region of human TBCCD1.		
Target	TBCCD1		
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.		
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		
Note	For research use only		
Application notes	For WB starting dilution is: 1:1000		
Isotype	Rabbit Ig		
Clonality	Polyclonal		
MW	64 kDa		
Uniprot ID	Q9NVR7		
NCBI	Q9NVR7		
Dilution Range	For WB starting dilution is: 1:1000		
Expiration Date	12 months from date of receipt.		