

Product Datasheet

LY6G6C Antibody (orb1270751)

Description

LY6G6C Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This LY6G6C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 62-90 amino acids from the C-terminal region of human LY6G6C.

Target

LY6G6C

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

14 kDa

Uniprot ID

[O95867](#)

NCBI

[O95867](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.

NCI-H292

72



Western blot analysis
in NCI-H292 cell l...



LY6G6C antibody
immunohistochemistry
ana...