

## Product Datasheet

### HTR3E Antibody (orb1270759)

**Description**

HTR3E Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

FC, WB

**Applications**
**Immunogen**

This 5HT3E antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 126-154 amino acids from the Central region of human 5HT3E.

**Target**

HTR3E

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**


51 kDa

**Uniprot ID**
**A5X5Y0**
**NCBI**
**A5X5Y0**
**Dilution Range**


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

**Expiration Date**

12 months from date of receipt.



Western blot analysis in CEM cell line l...



Flow cytometric analysis of CEM cells (r...