


Product Datasheet

MB21D1 Antibody (orb1270774)


Description

MB21D1 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	This CF150 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 266-295 amino acids from the Central region of human CF150.
Target	MB21D1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For FACS starting dilution is: 1:25For WB starting dilution is: 1:500-1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	Q8N884
NCBI	Q8N884
Dilution Range	For FACS starting dilution is: 1:25For WB starting dilution is: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Overlay histogram showing U-2OS cells st...



Flow cytometric analysis of A549 cells u...



Western blot analysis of lysate from THP...