

Product Datasheet

CYP2A7 Antibody (orb1270798)

Description

CYP2A7 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This CYP2A7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 121-148 amino acids from the Central region of human CYP2A7.

Target

CYP2A7

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


56 kDa

Uniprot ID
P20853
NCBI
P20853
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in MDA-MB231 cell ...



CYP2A7 antibody
immunohistochemistry
ana...



Flow cytometric
analysis of MDA-
MB231 ce...