

Product Datasheet

SMYD2 Antibody (orb1270820)

Description
SMYD2 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This SMYD2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 120-152 amino acids from the N-terminal region of human SMYD2.

Target

SMYD2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


50 kDa

Uniprot ID
[Q9NRG4](#)
NCBI
[Q9NRG4](#)
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Western blot
analysis of
SMYD2 using
rab...



Western blot
analysis in
NCI-H292 cell
l...