

Product Datasheet

KCNT2 Antibody (orb1270861)

Description
KCNT2 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This KCNT2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1069-1097 amino acids from the C-terminal region of human KCNT2.

Target

KCNT2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


131 kDa

Uniprot ID
Q6UVM3
NCBI
Q6UVM3
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.

MDA-MB231


Western blot analysis
in MDA-MB231 cell ...



KCNT2 antibody
immunohistochemistry
anal...



Flow cytometric
analysis of MDA-
MB231 ce...