

Product Datasheet

DNMT3A Antibody (orb1270864)

Description

DNMT3A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This Dnmt3a antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 457-486 amino acids from human Dnmt3a.

Target

DNMT3A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

102 kDa

Uniprot ID

[Q9Y6K1](#)

NCBI

[Q9Y6K1](#)

Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.

250

100

75

50

37


25

20


15

10

Western blot analysis of anti-Dnmt3a Pab...



Lysates from mice thymus tissue after ra...



Formalin-fixed and paraffin-embedded hum...