

Product Datasheet

ENTHD1 Antibody (orb1270867)

Description

ENTHD1 Antibody

130
95

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, WB

Applications

Immunogen

This ENTHD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 169-198 amino acids from the N-terminal region of human ENTHD1.

Target

ENTHD1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

68 kDa

Uniprot ID

[Q8IYW4](#)

NCBI


[Q8IYW4](#)

Dilution Range


For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis in A2058 cell line...



Flow cytometric analysis of A2058 cells ...