

Product Datasheet

ZDHHC1 Antibody (orb1270873)

Description

ZDHHC1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This ZDHHC1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 14-43 amino acids from the N-terminal region of human ZDHHC1.

Target

ZDHHC1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

[Q8WTX9](#)

NCBI

[Q8WTX9](#)


Dilution Range

For WB starting dilution is: 1:1000


Expiration Date

12 months from date of receipt.

m.Neuro-2a



Western blot analysis in mouse Neuro-2a ...



Western blot analysis in HepG2 cell line...