


Product Datasheet

Musculin antibody (orb127088)

Description	Rabbit Polyclonal to Musculin.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O60682
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...