

Product Datasheet

NNMT Antibody (orb1270895)

Description

NNMT Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This NNMT antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 77-106 amino acids from the Central region of human NNMT.

Target

NNMT

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

30 kDa

Uniprot ID

P40261

NCBI


P40261

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Western blot
analysis of
NNMT using
rabb...



Formalin-
fixed and
paraffin-
embedded
hum...