

Product Datasheet

RADIL Antibody (orb1270899)

Description

RADIL Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This RADIL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1044-1071 amino acids from the C-terminal region of human RADIL.

Target

RADIL

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


117 kDa

Uniprot ID
Q96JH8
NCBI
Q96JH8
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.

SK-BR-3
250


Western blot analysis
in SK-BR-3 cell li...



Western blot analysis
in mouse cerebellu...



RADIL antibody
immunohistochemistry
anal...