



Product Datasheet

alpha 2 Macroglobulin antibody (orb12709)

Description	Rabbit polyclonal antibody to alpha 2 Macroglobulin
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Alpha-2-Macroglobulin (Alpha-2-Macroglobulin (A2M))(NP_000005.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	163kDa
Uniprot ID	P01023
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of 293T cell...



Western blot analysis of extracts of var...