

Product Datasheet

ARL8A Antibody (orb1270936)

Description

ARL8A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This ARL8A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 116-142 amino acids from the Central region of human ARL8A.

Target

ARL8A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100
For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW




21 kDa

Uniprot ID
[Q96BM9](#)
NCBI
[Q96BM9](#)
Dilution Range

For WB starting dilution is: 1:1000
For IHC-P starting dilution is: 1:50~100
For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.


 Western blot analysis
in A375 cell line ...

 Western blot analysis
in mouse NIH-3T3 c...

 ARL8A Antibody
immunohistochemistry
anal...