

Product Datasheet

PPP1R14C Antibody (orb1270940)

Description

PPP1R14C Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

 Western blot analysis
in mouse heart tis...

Immunogen

This PPP1R14C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 126-155 amino acids from the C-terminal region of human PPP1R14C.

Target

PPP1R14C

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



18 kDa

Uniprot ID
[Q8TAE6](#)
NCBI
[Q8TAE6](#)
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date

12 months from date of receipt.


 PPP1R14C Antibody
immunohistochemistry
an...