

## Product Datasheet

### TFB2M Antibody (orb1270943)

## Description

TFB2M Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This TFB2M antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 368-396 amino acids from the C-terminal region of human TFB2M.

### Target

TFB2M

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

45 kDa

### Uniprot ID

[Q9H5Q4](#)

### NCBI

[Q9H5Q4](#)

### Dilution Range


For IHC-P starting dilution is: 1:50~100For WB starting dilution is: 1:1000

### Expiration Date

12 months from date of receipt.



TFB2M Antibody immunohistochemistry anal...



Western blot analysis in K562, WiDr cell...