

## Product Datasheet

### CYP3A4 Antibody (orb1271016)

## Description

CYP3A4 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested Applications

FC, IHC-P, WB

### Immunogen

This CYP3A4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 228-255 amino acids from the Central region of human CYP3A4.

### Target

CYP3A4

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

57 kDa

### Uniprot ID

[P08684](#)

### NCBI


[P08684](#)

### Dilution Range


For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

### Expiration Date


12 months from date of receipt.



Flow cytometric analysis of CEM cells (r...)



Western blot analysis of lysate from hum...



Western blot analysis of lysates from He...