

Product Datasheet

LGR5 Antibody (orb1271033)

Description

LGR5 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

This LGR5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 689-719 amino acids from human LGR5.

Target

LGR5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

100 kDa

Uniprot ID

[O75473](#)

NCBI


[O75473](#)

Dilution Range


For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:2000

Expiration Date


12 months from date of receipt.



Overlay histogram showing SH-SY5Y cells ...



Overlay histogram showing SH-SY5Y cells ...



Overlay histogram showing HepG2 cells st...