

Product Datasheet

LIN28B Antibody (orb1271039)

Description

LIN28B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

This LIN28B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 95-128 amino acids from the Central region of human LIN28B.

Target

LIN28B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

27 kDa

Uniprot ID

[Q6ZN17](#)

NCBI


[Q6ZN17](#)

Dilution Range


For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Antibody staining
LIN28B in
human
testis...



Antibody staining
LIN28B in
human
testis...



Western Blot
at 1:2000
dilution +
NCCIT ...