

Product Datasheet

HTRA1 Antibody (orb1271044)

Description

HTRA1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This HtrA1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 116-147 amino acids from the N-terminal region of human HtrA1.

Target

HTRA1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

51 kDa

Uniprot ID

[Q92743](#)

NCBI


[Q92743](#)

Dilution Range


For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot at
1:2000 dilution +
mouse ...



HtrA1 Antibody
immunohistochemistry
anal...



Confocal
immunofluorescent
analysis of H...