

Product Datasheet

Foxp3 Antibody (PE) (orb1271103)

Description	Foxp3 Antibody (PE)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC
Target	Foxp3
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b
Clonality	Monoclonal
Uniprot ID	Q99JB6
Expiration Date	12 months from date of receipt.