

Product Datasheet

OCT3/4 Antibody [EM92] (PE) (orb1271245)

Description	OCT3/4 Antibody [EM92] (PE)
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	FC
Target	POU5F1; Pou5f1
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P20263 , Q01860
Expiration Date	12 months from date of receipt.