

Product Datasheet

IL17a Antibody (FITC) (orb1271264)

Description	IL17a Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse, Rat
Conjugation	FITC
Tested Applications	FC
Target	IL17a
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	Q62386, Q61453
Expiration Date	12 months from date of receipt.