

Product Datasheet

TLR4 Antibody (PE) (orb1271301)

Description	TLR4 Antibody (PE)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Target	TLR4
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	O00206
Expiration Date	12 months from date of receipt.