

## Product Datasheet

### Cd40Ig Antibody (PE) (orb1271392)

|                            |                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Cd40Ig Antibody (PE)                                                                                                   |
| <b>Species/Host</b>        | Hamster                                                                                                                |
| <b>Reactivity</b>          | Mouse                                                                                                                  |
| <b>Conjugation</b>         | PE                                                                                                                     |
| <b>Tested Applications</b> | FC                                                                                                                     |
| <b>Target</b>              | Cd40Ig                                                                                                                 |
| <b>Preservatives</b>       | Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.              |
| <b>Form/Appearance</b>     | liquid                                                                                                                 |
| <b>Concentration</b>       | batch dependent                                                                                                        |
| <b>Storage</b>             | The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze. |
| <b>Note</b>                | For research use only                                                                                                  |
| <b>Isotype</b>             | Armenian Hamster IgG                                                                                                   |
| <b>Clonality</b>           | Monoclonal                                                                                                             |
| <b>Uniprot ID</b>          | <b>P27548</b>                                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                        |