

Product Datasheet

Itgb1 Antibody (PE) (orb1271461)

Description	Itgb1 Antibody (PE)
Species/Host	Hamster
Reactivity	Mouse, Rat
Conjugation	PE
Tested Applications	FC
Target	Itgb1
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Uniprot ID	P09055
Expiration Date	12 months from date of receipt.