


Product Datasheet


PRKAG2 antibody (orb127161)

Description	Rabbit Polyclonal to PRKAG2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UGJ0
Dilution Range	IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Peptide - +

Immunohistochemical analysis of formalin...



Peptide - +

Confocal immunofluorescence analysis of ...