

Product Datasheet

CD16/CD32 Antibody [2.4G2] (orb1271675)

Description	CD16/CD32 Antibody [2.4G2]
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC, IF, IP
Target	Fcgr3; Fcgr2b
Preservatives	Phosphate-buffered aqueous solution, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C . Do not freeze.
Note	For research use only
Isotype	Rat IgG2b
Clonality	Monoclonal
Uniprot ID	P08101, P08508
Expiration Date	12 months from date of receipt.