


Product Datasheet

Ptpcr Antibody (FITC) (orb1271696)

Description	Ptpcr Antibody (FITC)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Target	Ptpcr
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P06800
Expiration Date	12 months from date of receipt.



C57Bl/6
bone marrow
cells were
stained w...