


Product Datasheet

Ptprc Antibody (FITC) (orb1271711)

Description	Ptprc Antibody (FITC)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Target	Ptprc
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P06800
Expiration Date	12 months from date of receipt.



SJL
splenocytes
were stained
with FITC A...