

Product Datasheet

CD44 Antibody (BioWave Violet 450) (orb1271836)

Description	CD44 Antibody (BioWave Violet 450)
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	BioWave Violet 450
Tested Applications	FC
Target	CD44
Preservatives	Phosphate-buffered aqueous solution, $\leq 0.09\%$ Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P16070, P15379
Expiration Date	12 months from date of receipt.