

Product Datasheet

CD44 Antibody (orb1271840)

Description	CD44 Antibody
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC
Target	CD44
Preservatives	Phosphate-buffered aqueous solution, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C . Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P16070, P15379
Expiration Date	12 months from date of receipt.