



## Product Datasheet

### BIRC2 antibody (orb127187)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to BIRC2.               |
| <b>Species/Host</b>        | Rabbit                                             |
| <b>Reactivity</b>          | Human, Mouse                                       |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | IHC, WB                                            |
| <b>Immunogen</b>           | N-terminal                                         |
| <b>Concentration</b>       | 1mg/ml                                             |
| <b>Note</b>                | For research use only                              |
| <b>Clonality</b>           | Polyclonal                                         |
| <b>Uniprot ID</b>          | <b>P15927</b>                                      |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:1000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |



Western blot analysis of extracts from H...



Immunohistochemical analysis of formalin...