


Product Datasheet

Cd3e Antibody (PerCP/Cy5.5) (orb1271974)

Description	Cd3e Antibody (PerCP/Cy5.5)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Target	Cd3e
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Uniprot ID	P22646
Expiration Date	12 months from date of receipt.



C57Bl/6
splenocytes
were stained
with Pe...