



Product Datasheet

S6K-alpha2 antibody (orb127206)

Description	Rabbit Polyclonal to S6K-alpha2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15349
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from J...



Immunohistochemical analysis of formalin...