

Product Datasheet

Cd2 Antibody (FITC) (orb1272060)

Description	Cd2 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Target	Cd2
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, lambda
Clonality	Monoclonal
Uniprot ID	P08920
Expiration Date	12 months from date of receipt.